

Antibodies Product Datasheet

Product Name:	Influenza A virus Nucleoprotein Monoclonal Antibody
Catalog Number:	N0374
Reactivity:	A/Korea/426/68(H2N2)
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Influenza A virus nucleoprotein expressed in E. Coli strain BL21(DE3)
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Influenza A virus Nucleoprotein Monoclonal Antibody detects Influenza A virus Nucleoprotein protein.
Storage and Stability:	-20°C/1 year