

Antibodies Product Datasheet

Product Name:	Integrin α 5 Monoclonal Antibody
Catalog Number:	N0378
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Integrin α 5 expressed in E. Coli.
Concentration:	
Gene Name:	ITGA5
Protein Name:	Integrin alpha-5
Other Name:	ITGA5; FNRA; Integrin alpha-5; CD49 antigen-like family member E; Fibronectin receptor subunit alpha; Integrin alpha-F; VLA-5; CD antigen CD49e
Human Swiss Prot No.:	P08648
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Integrin α 5 Monoclonal Antibody detects endogenous levels of Integrin α 5 protein.
Storage and Stability:	-20°C/1 year