

Antibodies Product Datasheet

Product Name:	IRAK-4 Monoclonal Antibody
Catalog Number:	N0379
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF (paraffin section), FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human IRAK-4 expressed in E. Coli.
Concentration:	
Gene Name:	IRAK4
Protein Name:	Interleukin-1 receptor-associated kinase 4
Other Name:	
Human Swiss Prot No.:	Q8TDF7
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	IRAK-4 Monoclonal Antibody detects endogenous levels of IRAK-4 protein.
Storage and Stability:	-20°C/1 year