

Antibodies Product Datasheet

Product Name:	JAK2 Monoclonal Antibody
Catalog Number:	N0384
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of JAK2 (745-955aa) expressed in E. Coli.
Concentration:	
Gene Name:	JAK2
Protein Name:	Tyrosine-protein kinase JAK2
Other Name:	JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2
Human Swiss Prot No.:	O60674
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	JAK2 Monoclonal Antibody detects endogenous levels of JAK2 protein.
Storage and Stability:	-20°C/1 year