

Antibodies Product Datasheet

Product Name:	JAK3 Monoclonal Antibody
Catalog Number:	N0385
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human JAK3 expressed in E. Coli.
Concentration:	
Gene Name:	JAK3
Protein Name:	Tyrosine-protein kinase JAK3
Other Name:	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
Human Swiss Prot No.:	P52333
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	JAK3 Monoclonal Antibody detects endogenous levels of JAK3 protein.
Storage and Stability:	-20°C/1 year