

Antibodies Product Datasheet

Product Name:	JMJD1A Monoclonal Antibody
Catalog Number:	N0387
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human JMJD1A expressed in E. Coli.
Concentration:	
Gene Name:	KDM3A
Protein Name:	Lysine-specific demethylase 3A
Other Name:	KDM3A; JHDM2A; JMJD1; JMJD1A; KIAA0742; TSGA; Lysine-specific demethylase 3A; JmjC domain-containing histone demethylation protein 2A; Jumonji domain-containing protein 1A
Human Swiss Prot No.:	Q9Y4C1
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	JMJD1A Monoclonal Antibody detects endogenous levels of JMJD1A protein.
Storage and Stability:	-20°C/1 year