

Antibodies Product Datasheet

Product Name:	JMJD2A Monoclonal Antibody
Catalog Number:	N0388
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human JMJD2A expressed in E. Coli.
Concentration:	
Gene Name:	KDM4A
Protein Name:	Lysine-specific demethylase 4A
Other Name:	KDM4A; JHDM3A; JMJD2; JMJD2A; KIAA0677; Lysine-specific demethylase 4A; JmjC domain-containing histone demethylation protein 3A; Jumonji domain-containing protein 2A
Human Swiss Prot No.:	O75164
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	JMJD2A Monoclonal Antibody detects endogenous levels of JMJD2A protein.
Storage and Stability:	-20°C/1 year