

Antibodies Product Datasheet

Product Name:	JNK3 Monoclonal Antibody
Catalog Number:	N0389
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human JNK3 (aa28-233) expressed in E. Coli.
Concentration:	
Gene Name:	MAPK10
Protein Name:	Mitogen-activated protein kinase 10
Other Name:	MAPK10; JNK3; JNK3A; PRKM10; SAPK1B; Mitogen-activated protein kinase 10; MAP kinase 10; MAPK 10; MAP kinase p49 3F12; Stress-activated protein kinase 1b; SAPK1b; Stress-activated protein kinase JNK3; c-Jun N-terminal kinase 3
Human Swiss Prot No.:	P53779
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	JNK3 Monoclonal Antibody detects endogenous levels of JNK3 protein.
Storage and Stability:	-20°C/1 year