

Antibodies Product Datasheet

Product Name:	KID Monoclonal Antibody
Catalog Number:	N0394
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KID expressed in E. Coli.
Concentration:	
Gene Name:	KIF22
Protein Name:	Kinesin-like protein KIF22
Other Name:	KIF22; KID; KNSL4; Kinesin-like protein KIF22; Kinesin-like DNA-binding protein; Kinesin-like protein 4
Human Swiss Prot No.:	Q14807
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	KID Monoclonal Antibody detects endogenous levels of KID protein.
Storage and Stability:	-20°C/1 year