

Antibodies Product Datasheet

Product Name:	KLHL1 Monoclonal Antibody
Catalog Number:	N0396
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KLHL1 expressed in E. Coli.
Concentration:	
Gene Name:	KLHL1
Protein Name:	Kelch-like protein 1
Other Name:	KLHL1; KIAA1490; Kelch-like protein 1
Human Swiss Prot No.:	Q9NR64
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	KLHL1 Monoclonal Antibody detects endogenous levels of KLHL1 protein.
Storage and Stability:	-20°C/1 year