

Antibodies Product Datasheet

Product Name:	KLHL11 Monoclonal Antibody
Catalog Number:	N0397
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KLHL11 expressed in E. Coli.
Concentration:	
Gene Name:	KLHL11
Protein Name:	Kelch-like protein 11
Other Name:	KLHL11; Kelch-like protein 11
Human Swiss Prot No.:	Q9NVR0
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	KLHL11 Monoclonal Antibody detects endogenous levels of KLHL11 protein.
Storage and Stability:	-20°C/1 year