

## Antibodies Product Datasheet

|                               |                                                                                       |
|-------------------------------|---------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | KLHL12 Monoclonal Antibody                                                            |
| <b>Catalog Number:</b>        | N0399                                                                                 |
| <b>Reactivity:</b>            | Human                                                                                 |
| <b>Applications:</b>          | WB,ELISA                                                                              |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.   |
| <b>Modification:</b>          | Unmodified/ Total                                                                     |
| <b>Source:</b>                | Monoclonal Mouse                                                                      |
| <b>Purification:</b>          | Affinity purification                                                                 |
| <b>Immunogen:</b>             | Purified recombinant fragment of human KLHL12 expressed in E. Coli.                   |
| <b>Concentration:</b>         |                                                                                       |
| <b>Gene Name:</b>             | KLHL12                                                                                |
| <b>Protein Name:</b>          | Kelch-like protein 12                                                                 |
| <b>Other Name:</b>            | KLHL12; C3IP1; Kelch-like protein 12; CUL3-interacting protein 1; DKIR homolog; hDKIR |
| <b>Human Swiss Prot No.:</b>  | Q53G59                                                                                |
| <b>Formulation:</b>           | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.                    |
| <b>Specificity:</b>           | KLHL12 Monoclonal Antibody detects endogenous levels of KLHL12 protein.               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                          |