

Antibodies Product Datasheet

Product Name:	KLHL13 Monoclonal Antibody
Catalog Number:	N0400
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KLHL13 expressed in E. Coli.
Concentration:	
Gene Name:	KLHL13
Protein Name:	Kelch-like protein 13
Other Name:	KLHL13; BKLHD2; KIAA1309; Kelch-like protein 13; BTB and kelch domain-containing protein 2
Human Swiss Prot No.:	Q9P2N7
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	KLHL13 Monoclonal Antibody detects endogenous levels of KLHL13 protein.
Storage and Stability:	-20°C/1 year