

Antibodies Product Datasheet

Product Name:	KLHL25 Monoclonal Antibody
Catalog Number:	N0403
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KLHL25 expressed in E. Coli.
Concentration:	
Gene Name:	KLHL25
Protein Name:	Ectoderm-neural cortex protein 2
Other Name:	KLHL25; ENC2; Ectoderm-neural cortex protein 2; ENC-2; Kelch-like protein 25
Human Swiss Prot No.:	Q9H0H3
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	KLHL25 Monoclonal Antibody detects endogenous levels of KLHL25 protein.
Storage and Stability:	-20°C/1 year