

Antibodies Product Datasheet

Product Name:	KSHV ORF62 Monoclonal Antibody
Catalog Number:	N0407
Reactivity:	KsheepV
Applications:	WB, ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunocytochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human KSHV ORF62 expressed in E. Coli.
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	KSHV ORF62 Monoclonal Antibody detects KSHV ORF62 protein.
Storage and Stability:	-20°C/1 year