

Antibodies Product Datasheet

Product Name:	LAL Monoclonal Antibody
Catalog Number:	N0409
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of LAL expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	LIPA
Protein Name:	Lysosomal acid lipase/cholesteryl ester hydrolase
Other Name:	LIPA; Lysosomal acid lipase/cholesteryl ester hydrolase; Acid cholesteryl ester hydrolase; LAL; Cholesteryl esterase; Lipase A; Sterol esterase
Human Swiss Prot No.:	P38571
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	LAL Monoclonal Antibody detects endogenous levels of LAL protein.
Storage and Stability:	-20°C/1 year