

Antibodies Product Datasheet

Product Name:	Laminin β -1 Monoclonal Antibody
Catalog Number:	N0411
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Laminin β -1 (aa31-270) expressed in E. Coli.
Concentration:	
Gene Name:	LAMB1
Protein Name:	Laminin subunit beta-1
Other Name:	LAMB1; Laminin subunit beta-1; Laminin B1 chain; Laminin-1 subunit beta; Laminin-10 subunit beta; Laminin-12 subunit beta; Laminin-2 subunit beta; Laminin-6 subunit beta; Laminin-8 subunit beta
Human Swiss Prot No.:	P07942
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Laminin β -1 Monoclonal Antibody detects endogenous levels of Laminin β -1 protein.
Storage and Stability:	-20°C/1 year