

Antibodies Product Datasheet

Product Name:	Lck Monoclonal Antibody
Catalog Number:	N0412
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Lck expressed in E. Coli.
Concentration:	
Gene Name:	LCK
Protein Name:	Tyrosine-protein kinase Lck
Other Name:	LCK; Tyrosine-protein kinase Lck; Leukocyte C-terminal Src kinase; LSK; Lymphocyte cell-specific protein-tyrosine kinase; Protein YT16; Proto-oncogene Lck; T cell-specific protein-tyrosine kinase; p56-LCK
Human Swiss Prot No.:	P06239
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Lck Monoclonal Antibody detects endogenous levels of Lck protein.
Storage and Stability:	-20°C/1 year