

Antibodies Product Datasheet

Product Name:	LEDGF Monoclonal Antibody
Catalog Number:	N0413
Reactivity:	Human,Rat,Monkey
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human LEDGF expressed in E. Coli.
Concentration:	
Gene Name:	PSIP1
Protein Name:	PC4 and SFRS1-interacting protein
Other Name:	PSIP1; DFS70; LEDGF; PSIP2; PC4 and SFRS1-interacting protein; CLL-associated antigen KW-7; Dense fine speckles 70 kDa protein; DFS 70; Lens epithelium-derived growth factor; Transcriptional coactivator p75/p52
Human Swiss Prot No.:	O75475
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	LEDGF Monoclonal Antibody detects endogenous levels of LEDGF protein.
Storage and Stability:	-20°C/1 year