

Antibodies Product Datasheet

Product Name:	LHX2 Monoclonal Antibody
Catalog Number:	N0414
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human LHX2 expressed in E. Coli.
Concentration:	
Gene Name:	LHX2
Protein Name:	LIM/homeobox protein Lhx2
Other Name:	LHX2; LH2; LIM/homeobox protein Lhx2; Homeobox protein LH-2; LIM homeobox protein 2
Human Swiss Prot No.:	P50458
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	LHX2 Monoclonal Antibody detects endogenous levels of LHX2 protein.
Storage and Stability:	-20°C/1 year