

Antibodies Product Datasheet

Product Name:	LIN-28 Monoclonal Antibody
Catalog Number:	N0415
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of LIN-28 (aa93-209) expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	LIN28A
Protein Name:	Protein lin-28 homolog A
Other Name:	LIN28A; CSDD1; LIN28; ZCCHC1; Protein lin-28 homolog A; Lin-28A; Zinc finger CCHC domain-containing protein 1
Human Swiss Prot No.:	Q9H9Z2
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in 0.01M Phosphate buffer, pH 7.4 containing 0.03% sodium azide.
Specificity:	LIN-28 Monoclonal Antibody detects endogenous levels of LIN-28 protein.
Storage and Stability:	-20°C/1 year