

Antibodies Product Datasheet

Product Name:	LMP7A Monoclonal Antibody
Catalog Number:	N0418
Reactivity:	Human,Rat
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human LMP7A expressed in E. Coli.
Concentration:	
Gene Name:	PSMB8
Protein Name:	Proteasome subunit beta type-8
Other Name:	PSMB8; LMP7; PSMB5i; RING10; Y2; Proteasome subunit beta type-8; Low molecular mass protein 7; Macropain subunit C13; Multicatalytic endopeptidase complex subunit C13; Proteasome component C13; Proteasome subunit beta-5i; Really interesting
Human Swiss Prot No.:	P28062
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	LMP7A Monoclonal Antibody detects endogenous levels of LMP7A protein.
Storage and Stability:	-20°C/1 year