

Antibodies Product Datasheet

Product Name:	Lsk Monoclonal Antibody
Catalog Number:	N0422
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Lsk expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	MATK
Protein Name:	Megakaryocyte-associated tyrosine-protein kinase
Other Name:	MATK; CTK; HYL; Megakaryocyte-associated tyrosine-protein kinase; CSK homologous kinase; CHK; Hematopoietic consensus tyrosine-lacking kinase; Protein kinase HYL; Tyrosine-protein kinase CTK
Human Swiss Prot No.:	P42679
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 50% glycerol and 0.03% sodium azide.
Specificity:	Lsk Monoclonal Antibody detects endogenous levels of Lsk protein.
Storage and Stability:	-20°C/1 year