

## Antibodies Product Datasheet

|                               |                                                                                                                                                                 |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Lyn Monoclonal Antibody                                                                                                                                         |
| <b>Catalog Number:</b>        | N0423                                                                                                                                                           |
| <b>Reactivity:</b>            | Human                                                                                                                                                           |
| <b>Applications:</b>          | WB,ELISA                                                                                                                                                        |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.                                                                             |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                               |
| <b>Source:</b>                | Monoclonal Mouse                                                                                                                                                |
| <b>Purification:</b>          | Affinity purification                                                                                                                                           |
| <b>Immunogen:</b>             | Purified recombinant fragment of Lyn expressed in E. Coli.                                                                                                      |
| <b>Concentration:</b>         |                                                                                                                                                                 |
| <b>Gene Name:</b>             | LYN                                                                                                                                                             |
| <b>Protein Name:</b>          | Tyrosine-protein kinase Lyn                                                                                                                                     |
| <b>Other Name:</b>            | LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn |
| <b>Human Swiss Prot No.:</b>  | P07948                                                                                                                                                          |
| <b>Formulation:</b>           | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.                                                                                              |
| <b>Specificity:</b>           | Lyn Monoclonal Antibody detects endogenous levels of Lyn protein.                                                                                               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                    |