

Antibodies Product Datasheet

Product Name:	Mcl-1 Monoclonal Antibody
Catalog Number:	N0429
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human MCL-1 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	MCL1
Protein Name:	Induced myeloid leukemia cell differentiation protein Mcl-1
Other Name:	MCL1; BCL2L3; Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT
Human Swiss Prot No.:	Q07820
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	Mcl-1 Monoclonal Antibody detects endogenous levels of Mcl-1 protein.
Storage and Stability:	-20°C/1 year