

Antibodies Product Datasheet

Product Name:	MCP-1 Monoclonal Antibody
Catalog Number:	N0430
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human MCP-1 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	CCL2
Protein Name:	C-C motif chemokine 2
Other Name:	CCL2; MCP1; SCYA2; C-C motif chemokine 2; HC11; Monocyte chemoattractant protein 1; Monocyte chemotactic and activating factor; MCAF; Monocyte chemotactic protein 1; MCP-1; Monocyte secretory protein JE; Small-inducible cytokine A2
Human Swiss Prot No.:	P13500
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	MCP-1 Monoclonal Antibody detects endogenous levels of MCP-1 protein.
Storage and Stability:	-20°C/1 year