

Antibodies Product Datasheet

Product Name:	M-CSF Monoclonal Antibody
Catalog Number:	N0431
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human M-CSF expressed in E. Coli.
Concentration:	
Gene Name:	CSF1
Protein Name:	Macrophage colony-stimulating factor 1
Other Name:	CSF1; Macrophage colony-stimulating factor 1; CSF-1; M-CSF; MCSF; Lanimostim
Human Swiss Prot No.:	P09603
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	M-CSF Monoclonal Antibody detects endogenous levels of M-CSF protein.
Storage and Stability:	-20°C/1 year