

Antibodies Product Datasheet

Product Name:	MEK Kinase-2 Monoclonal Antibody
Catalog Number:	N0433
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human MEK Kinase-2 expressed in E. Coli.
Concentration:	
Gene Name:	MAP3K2
Protein Name:	Mitogen-activated protein kinase kinase kinase 2
Other Name:	MAP3K2; MAPKKK2; MEKK2; Mitogen-activated protein kinase kinase kinase 2; MAPK/ERK kinase kinase 2; MEK kinase 2; MEKK 2
Human Swiss Prot No.:	Q9Y2U5
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MEK Kinase-2 Monoclonal Antibody detects endogenous levels of MEK Kinase-2 protein.
Storage and Stability:	-20°C/1 year