

Antibodies Product Datasheet

Product Name:	MerTK Monoclonal Antibody
Catalog Number:	N0439
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of MerTK expressed in E. Coli.
Concentration:	
Gene Name:	MERTK
Protein Name:	Tyrosine-protein kinase Mer
Other Name:	MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK
Human Swiss Prot No.:	Q12866
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MerTK Monoclonal Antibody detects endogenous levels of MerTK protein.
Storage and Stability:	-20°C/1 year