

Antibodies Product Datasheet

Product Name:	Mi2- α Monoclonal Antibody
Catalog Number:	N0441
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Mi2- α expressed in E. Coli.
Concentration:	
Gene Name:	CHD3
Protein Name:	Chromodomain-helicase-DNA-binding protein 3
Other Name:	CHD3; Chromodomain-helicase-DNA-binding protein 3; CHD-3; ATP-dependent helicase CHD3; Mi-2 autoantigen 240 kDa protein; Mi2-alpha; Zinc finger helicase; hZFH
Human Swiss Prot No.:	Q12873
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Mi2- α Monoclonal Antibody detects endogenous levels of Mi2- α protein.
Storage and Stability:	-20°C/1 year