

Antibodies Product Datasheet

Product Name:	MOF Monoclonal Antibody
Catalog Number:	N0447
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human MOF expressed in E. Coli.
Concentration:	
Gene Name:	KAT8
Protein Name:	Histone acetyltransferase KAT8
Other Name:	KAT8; MOF; MYST1; PP7073; Histone acetyltransferase KAT8; Lysine acetyltransferase 8; MOZ; YBF2/SAS3, SAS2 and TIP60 protein 1; MYST-1; hMOF
Human Swiss Prot No.:	Q9H7Z6
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MOF Monoclonal Antibody detects endogenous levels of MOF protein.
Storage and Stability:	-20°C/1 year