

Antibodies Product Datasheet

Product Name:	Msi1 Monoclonal Antibody
Catalog Number:	N0450
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Msi1 expressed in E. Coli.
Concentration:	
Gene Name:	MSI1
Protein Name:	RNA-binding protein Musashi homolog 1
Other Name:	MSI1; RNA-binding protein Musashi homolog 1; Musashi-1
Human Swiss Prot No.:	O43347
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Msi1 Monoclonal Antibody detects endogenous levels of Msi1 protein.
Storage and Stability:	-20°C/1 year