

Antibodies Product Datasheet

Product Name:	Msx-1 Monoclonal Antibody
Catalog Number:	N0452
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Msx-1 expressed in E. Coli.
Concentration:	
Gene Name:	MSX1
Protein Name:	Homeobox protein MSX-1
Other Name:	MSX1; HOX7; Homeobox protein MSX-1; Homeobox protein Hox-7; Msh homeobox 1-like protein
Human Swiss Prot No.:	P28360
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Msx-1 Monoclonal Antibody detects endogenous levels of Msx-1 protein.
Storage and Stability:	-20°C/1 year