

Antibodies Product Datasheet

Product Name:	MTDH Monoclonal Antibody
Catalog Number:	N0453
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human MTDH expressed in E. Coli.
Concentration:	
Gene Name:	MTDH
Protein Name:	Protein LYRIC
Other Name:	MTDH; AEG1; LYRIC; Protein LYRIC; 3D3/LYRIC; Astrocyte elevated gene-1 protein; AEG-1; Lysine-rich CEACAM1 co-isolated protein; Metadherin; Metastasis adhesion protein
Human Swiss Prot No.:	Q86UE4
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MTDH Monoclonal Antibody detects endogenous levels of MTDH protein.
Storage and Stability:	-20°C/1 year