

Antibodies Product Datasheet

Product Name:	MuSK Monoclonal Antibody
Catalog Number:	N0456
Reactivity:	Human
Applications:	IHC-p,IF/ICC,ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human MuSK (aa24-209) fused with hIgGFc tag expressed in HEK293 cell line.
Concentration:	
Gene Name:	MUSK
Protein Name:	Muscle, skeletal receptor tyrosine-protein kinase
Other Name:	MUSK; Muscle; skeletal receptor tyrosine-protein kinase; Muscle-specific tyrosine-protein kinase receptor; MuSK; Muscle-specific kinase receptor
Human Swiss Prot No.:	O15146
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MuSK Monoclonal Antibody detects endogenous levels of MuSK protein.
Storage and Stability:	-20°C/1 year