

Antibodies Product Datasheet

Product Name:	Myoglobin Monoclonal Antibody
Catalog Number:	N0461
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Myoglobin expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	MB
Protein Name:	Myoglobin
Other Name:	MB; Myoglobin
Human Swiss Prot No.:	P02144
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	Myoglobin Monoclonal Antibody detects endogenous levels of Myoglobin protein.
Storage and Stability:	-20°C/1 year