

Antibodies Product Datasheet

Product Name:	NAC-1 Monoclonal Antibody
Catalog Number:	N0462
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human NAC-1 expressed in E. Coli.
Concentration:	
Gene Name:	BTB14
Protein Name:	Nucleus accumbens-associated protein
Other Name:	NACC1; BTBD14B; NAC1; Nucleus accumbens-associated protein 1; NAC-1; BTB/POZ domain-containing protein 14B
Human Swiss Prot No.:	Q96RE7
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	NAC-1 Monoclonal Antibody detects endogenous levels of NAC-1 protein.
Storage and Stability:	-20°C/1 year