

Antibodies Product Datasheet

Product Name:	N-cadherin Monoclonal Antibody
Catalog Number:	N0464
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human N-cadherin expressed in E. Coli.
Concentration:	
Gene Name:	CDH2
Protein Name:	Cadherin-2
Other Name:	CDH2; CDHN; NCAD; Cadherin-2; CDw325; Neural cadherin; N-cadherin; CD antigen CD325
Human Swiss Prot No.:	P19022
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	N-cadherin Monoclonal Antibody detects endogenous levels of N-cadherin protein.
Storage and Stability:	-20°C/1 year