

Antibodies Product Datasheet

Product Name:	N-CoR Monoclonal Antibody
Catalog Number:	N0466
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of N-CoR (aa1-192) expressed in E. Coli.
Concentration:	
Gene Name:	NCOR1
Protein Name:	Nuclear receptor corepressor 1
Other Name:	NCOR1; KIAA1047; Nuclear receptor corepressor 1; N-CoR; N-CoR1
Human Swiss Prot No.:	O75376
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	N-CoR Monoclonal Antibody detects endogenous levels of N-CoR protein.
Storage and Stability:	-20°C/1 year