

Antibodies Product Datasheet

Product Name:	Neurogenin 3 Monoclonal Antibody
Catalog Number:	N0470
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Neurogenin-3 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	NEUROG3
Protein Name:	Neurogenin-3
Other Name:	NEUROG3; ATOH5; BHLHA7; NGN3; Neurogenin-3; NGN-3; Class A basic helix-loop-helix protein 7; bHLHa7; Protein atonal homolog 5
Human Swiss Prot No.:	Q9Y4Z2
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 50% glycerol and 0.03% sodium azide.
Specificity:	Neurogenin 3 Monoclonal Antibody detects endogenous levels of Neurogenin 3 protein.
Storage and Stability:	-20°C/1 year

N0470 | Neurogenin 3 Monoclonal Antibody

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022