

Antibodies Product Datasheet

Product Name:	Nkx-2.5 Monoclonal Antibody
Catalog Number:	N0475
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Nkx-2.5 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	NKX2-5
Protein Name:	Homeobox protein Nkx-2.5
Other Name:	NKX2-5; CSX; NKX2.5; NKX2E; Homeobox protein Nkx-2.5; Cardiac-specific homeobox; Homeobox protein CSX; Homeobox protein NK-2 homolog E
Human Swiss Prot No.:	P52952
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 0.03% sodium azide.
Specificity:	Nkx-2.5 Monoclonal Antibody detects endogenous levels of Nkx-2.5 protein.
Storage and Stability:	-20°C/1 year