

Antibodies Product Datasheet

Product Name:	NOK Monoclonal Antibody
Catalog Number:	N0478
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of NOK expressed in E. Coli.
Concentration:	
Gene Name:	STYK1
Protein Name:	Tyrosine-protein kinase STYK1
Other Name:	STYK1; NOK; Tyrosine-protein kinase STYK1; Novel oncogene with kinase domain; Protein PK-unique; Serine/threonine/tyrosine kinase 1
Human Swiss Prot No.:	Q6J9G0
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	NOK Monoclonal Antibody detects endogenous levels of NOK protein.
Storage and Stability:	-20°C/1 year