

Antibodies Product Datasheet

Product Name:	NOS3 Monoclonal Antibody
Catalog Number:	N0479
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human NOS3 expressed in E. Coli.
Concentration:	
Gene Name:	NOS3
Protein Name:	Nitric oxide synthase, endothelial
Other Name:	NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII
Human Swiss Prot No.:	P29474
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	NOS3 Monoclonal Antibody detects endogenous levels of NOS3 protein.
Storage and Stability:	-20°C/1 year