

Antibodies Product Datasheet

Product Name:	OCT1 Monoclonal Antibody
Catalog Number:	N0481
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human OCT1 expressed in E. Coli.
Concentration:	
Gene Name:	SLC22A1
Protein Name:	Solute carrier family 22 member 1
Other Name:	SLC22A1; OCT1; Solute carrier family 22 member 1; Organic cation transporter 1; hOCT1
Human Swiss Prot No.:	O15245
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	OCT1 Monoclonal Antibody detects endogenous levels of 41183 protein.
Storage and Stability:	-20°C/1 year