

Antibodies Product Datasheet

Product Name:	OPG Monoclonal Antibody
Catalog Number:	N0487
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human OPG expressed in E. Coli.
Concentration:	
Gene Name:	TNFRSF11B
Protein Name:	Tumor necrosis factor receptor superfamily member 11B
Other Name:	TNFRSF11B; OCIF; OPG; Tumor necrosis factor receptor superfamily member 11B; Osteoclastogenesis inhibitory factor; Osteoprotegerin
Human Swiss Prot No.:	O00300
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	OPG Monoclonal Antibody detects endogenous levels of OPG protein.
Storage and Stability:	-20°C/1 year