

Antibodies Product Datasheet

Product Name:	p16 Monoclonal Antibody
Catalog Number:	N0493
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF (paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of p16 expressed in E. Coli.
Concentration:	
Gene Name:	CDKN2A
Protein Name:	Cyclin-dependent kinase inhibitor 2A isoform 4
Other Name:	CDKN2A; CDKN2; MLM; Cyclin-dependent kinase inhibitor 2A; isoform 4; p14ARF; p19ARF
Human Swiss Prot No.:	Q8N726
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	p16 Monoclonal Antibody detects endogenous levels of p16 protein.
Storage and Stability:	-20°C/1 year