

Antibodies Product Datasheet

Product Name:	p27 Monoclonal Antibody
Catalog Number:	N0495
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human p27 expressed in E. Coli.
Concentration:	
Gene Name:	CDKN1B
Protein Name:	Cyclin-dependent kinase inhibitor 1B
Other Name:	CDKN1B; KIP1; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; p27Kip1
Human Swiss Prot No.:	P46527
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	p27 Monoclonal Antibody detects endogenous levels of p27 protein.
Storage and Stability:	-20°C/1 year