

Antibodies Product Datasheet

Product Name:	p38 β Monoclonal Antibody
Catalog Number:	N0497
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of p38 β (aa251-363) expressed in E. Coli.
Concentration:	
Gene Name:	MAPK11
Protein Name:	Mitogen-activated protein kinase 11
Other Name:	MAPK11; PRKM11; SAPK2; SAPK2B; Mitogen-activated protein kinase 11; MAP kinase 11; MAPK 11; Mitogen-activated protein kinase p38 beta; MAP kinase p38 beta; p38b; Stress-activated protein kinase 2b; SAPK2b; p38-2
Human Swiss Prot No.:	Q15759
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	p38 β Monoclonal Antibody detects endogenous levels of p38 β protein.
Storage and Stability:	-20°C/1 year