

Antibodies Product Datasheet

Product Name:	PAKy Monoclonal Antibody
Catalog Number:	N0501
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PAKy expressed in E. Coli.
Concentration:	
Gene Name:	PAK2
Protein Name:	Serine/threonine-protein kinase PAK 2
Other Name:	PAK2; Serine/threonine-protein kinase PAK 2; Gamma-PAK; PAK65; S6/H4 kinase; p21-activated kinase 2; PAK-2; p58
Human Swiss Prot No.:	Q13177
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	PAKy Monoclonal Antibody detects endogenous levels of PAKy protein.
Storage and Stability:	-20°C/1 year