

Antibodies Product Datasheet

Product Name:	PEA3 Monoclonal Antibody
Catalog Number:	N0513
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PEA3 (aa50-109) expressed in E. Coli.
Concentration:	
Gene Name:	ETV4
Protein Name:	ETS translocation variant 4
Other Name:	ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
Human Swiss Prot No.:	P43268
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PEA3 Monoclonal Antibody detects endogenous levels of PEA3 protein.
Storage and Stability:	-20°C/1 year